




EQUINOX

G-NAIL QUICK TIPS

G-NAIL QUICK TIPS – overall kit layout


1. Nail selection




2. Nail Assembly & Insertion




3. Femoral Neck Screw Insertion



4. Distal Locking Screw Insertion



G-NAIL QUICK TIPS




1. G-NAIL: NAIL SELECTION – Implant tray (tray 5) required

Steps 1. Surgeon will view X-ray to determine nail length and diameter.

Diameter Guide*:


- Female: 9mm, 10mm
- Male: 11mm, 12mm

*approximation only.




2. G-NAIL: PREPARATION – Tray 1 & 5 required

Used: (1) Aiming arm (Static)
(2) Insertion handle
(3) G-Nail (Tray 5)
(4) Locking bolt
(5) Disco Ball wrench
(9) Guide Bar (optional)




- Steps:
1. Connect the Aiming Arm **No. 1.** and Insertion Handle **No. 2.** Guide Bar **No. 9** may be attached to the top of the insertion handle at this time, if required – especially for tough bone.
 2. Align the Aiming Arm tooth into the G-Nail notch.
 3. Use the Locking Bolt **No. 4** to secure the top of the nail to the aiming arm, tightening with the Disco Ball Wrench **No. 5**

2. G-NAIL: PREPARATION



3. G-NAIL: DRILLING & NAIL INSERTION – Tray 1 required

Used: (6) 3.2mm Guide Wire
(7) Wire Guide (optional)
(8) Guide Sleeve
(9) 17 Cannulated Bone Drill
* Chuck Attachment Required




- Steps:**
1. The 3.2mm Guide Wire **No. 6** is inserted through the top of the greater trochanter (Wire Guide optional)
 2. 17mm Cannulated Bone Drill **No. 9** is inserted over the guide wire, through the Drill Guide **No. 8** and used to drill into the cancellous bone.

4. G-NAIL: INSERTION – Tray 1 required

Used: (9) Guide Bar


(37) Mallet (Hospitals should have a single wrapped mallet, if not, use mallet from Tray 4)



- Steps:
1. Insert the G-Nail into the pre-drilled space
 2. Guide Bar **No. 9** and Mallet **No. 37** may be needed to penetrate tough bone.
 3. Guide bar **No. 9** removed if used.


5. LAG SCREW: PREP – Tray 1 & 2 required

Used: (11) Buttress Nut
(12) Threaded Sleeve
(13) Yellow Guide Sleeve
(14) Yellow Trocar
(6) 3.2mm Guide Wire (Tray 1)
(16) Cannulated Depth Gauge



- Steps:
1. Screw the Buttress Nut **No. 11** onto the Threaded Sleeve **No. 12**, progress approximately 2/3 of the way down (**reverse thread**)
 2. Insert Yellow Trocar **No. 14** into Yellow Guide Sleeve **No. 13** and then the combined construct into the Threaded Sleeve **No. 12**. The construct is then placed through the femoral neck aiming guide of the Aiming Arm **No. 1**.
 3. The Mallet impacts the construct to puncture the lateral cortex of the proximal femur before the Yellow Trocar **No. 14** is removed.
 4. The 3.2mm guide wire **No. 6** is inserted 10mm from the articular surface and then the yellow Guide Sleeve **No. 13** is removed and the Threaded Sleeve **No. 12** progressed onto the lateral cortex wall.
 5. The Cannulated Depth Gauge **No. 16** then determines the length of lag screw required with the yellow ringed end placed over the end of the threaded sleeve.

5. LAG SCREW: PREP – Tray 1 & 2 required




6. LAG SCREW: DRILLING – Tray 2 required

Used: (17) Depth Stop
(18) 10.5mm Reamer
(19) 10.5mm Bone Drill

Patient side → Desired depth should be readable



Surgeon side





- Steps:**
1. Depth Stop **No. 17** is calibrated according to the measured length (length shown on the **PATIENT SIDE** – 85mm in this example).
 2. 10.5mm Bone Drill **No. 19** is progressed over the guide wire to penetrate the cortical bone of the lateral cortex, **if required**, and then removed.
 3. 10.5mm Reamer **No. 18** is progressed over guide wire to drill for Lag Screw insertion.

7.1. LAG SCREW: INSERTION – Tray 2 required

Used: (20) T Handle Screwdriver
(21) Capture Rod
(22) Pressuriser (optional)
(23) Universal Hex Screwdriver


* See the next page for a detailed schematic




- Steps:**
1. Insert Capture Rod **No. 21** into T-handle **No. 20**.
 2. Screw Capture Rod into Lag Screw to create stable construct and hand to Surgeon.
 3. Pressuriser **No. 22** (Optional when interfragmentary compression of the fracture fragments is required. The Pressuriser is screwed onto the T-Handle shaft and screwed clockwise).
 4. Ensure the green line on the T Handle Screwdriver is facing **cranially** and in line with the **green markings** when fully inserted.
 5. Universal Screwdriver **No. 23** to engage the internal locking mechanism that locks the Lag Screw to the G-Nail. The Surgeon may perform one quarter turn back to generate micromotion of the lag screw.

7.2. LAG SCREW: INSERTION – Tray 2 required

Ensure the green line on the T Handle Screwdriver is facing cranially and in line with the green marking on the threaded sleeve when fully inserted.




- The **internal locking mechanism** comes pre-screwed into the G-Nail.
- Use the Universal Screwdriver to engage the internal locking mechanism and lock the Lag Screw to the G-Nail.



Excessive sterilization of the threaded sleeve and T-Handle Screwdriver may remove the green markings. Take care to ensure the green marking (or groove) is always aligned with the marking/groove on the threaded sleeve and is facing cranially.

8. DISTAL LOCKING SCREW: DRILLING – Tray 3 required


Used: (24) Protection Sleeve
(25) Green Drill Sleeve
(26) Green Trocar
(27) 4.2mm Bone drill
(calibrated depth gauge)
(28) Depth gauge
(optional)



- Steps:**
1. Insert Green Trocar **No. 26** into Green Drill Sleeve **No. 25** and then the combined construct into the Protection Sleeve **No. 24** and thread through the aiming arm.
 2. The Green Trocar **No. 26** is then impacted with the mallet to puncture the lateral cortex of the femur and is then removed with the Green Drill Sleeve **No. 25**.
 3. 4.2mm Bone Drill **No. 27** measures the length of the Distal Locking Screw. The drill is QC, so no chuck is required.
 4. Depth Gauge **No. 28** (optional).

9. DISTAL LOCKING SCREW: INSERTION – Tray 3 & 5 required

Used: (29) 4.0mm Green Cannulated Hex Screwdriver
(30) Capture rod
(31) Distal locking screw (Tray 5)





- Steps:**
1. Capture Rod **No. 30** is inserted into the 4.0mm Green Cannulated Hex Screwdriver **No. 29**.
 2. Attach Capture Rod **No. 30** firmly to the head of the Distal Locking Screw.
 3. Hand to the surgeon who will thread this through the Aiming Arm and secure.

10. NAIL REMOVAL – Auxillary & Implant Removal Tray required

Used:


- (29) 4.0mm Green Cannulated Hex Screwdriver
- (30) Capture rod
- (23) Universal Hex Screwdriver
- (20) T Handle Screwdriver
- (21) Capture Rod
- (38) Nail Extraction Sleeve
- (39) Guide Rod
- (37) Mallet



- Steps:**
1. **Distal Locking Screw Removal:** Capture Rod **No. 30** is inserted into the 4.0mm Green Cannulated Hex Screwdriver **No. 29** and the locking screw is extracted.
 2. **Lag Screw Removal:** The Universal 4.5 Hex Screwdriver **No. 23** is used to unscrew the set screw first. Capture Rod **No. 21** is inserted into the T Handle Screwdriver **No. 20** and the lag screw is extracted.
 3. **Nail Removal:** The Nail Extraction Sleeve **No. 38** is attached to the proximal end of the nail and a Guide Rod **No. 39** is then attached to the Nail Extraction sleeve. A Mallet **No. 37** is then used to extract the nail. (attachments shown in the next page)

*Please note that **No. 29, 30, 23, 20 & 21** are all left loose on the Auxillary Tray and do not have a proper laser engraved spot.*

10. NAIL REMOVAL – Auxillary & Implant Removal Tray required



Explanation – *ABCD* Markings on the Buttress Nut



The four letters, A, B, C and D serve the following function:

- Indicates the correct direction of turning
 - Turning from A to D, advances the trocar onto the proximal cortex (clockwise)
- Provides a measure of the amount of advancement
 - One full revolution (360°) = 2.0mm of compression
 - 90° turn eg. from A to B = 0.5mm of compression
 - 1800 turn eg. A to C = 1.0mm of compression.