

WEGO

Blood Collection Tubes



EQUINOX

OUR COMMITMENT

The Weigao Group was established in 1988, and has since become a top medical supplies manufacturer in China with a workforce of over 10000; principally involved in the research and development ("R&D"), production and sales of branded disposable medical devices.

The Vacuum blood collection tube is one of Weigao Group's most experienced production lines, which has been manufactured since its foundation in 1992. Weigao Group has since introduced certified vacuum blood collection tube, production and quality control to market complete with CE marks, ISO9001, ISO13485 and GMP certificates.

All manufacturing processes are quality controlled and maintained to the highest standards; Weigao's production process includes the use of imported Canadian Husky Injection Mould machine. Our vacuum blood collection tube has included all additives and sharp blood collection needles, holder.

SERUM	PLASMA	COAGULATION	SERUM SEPARATOR	ESR	WHOLE BLOOD	BLOOD GLUCOSE
PLAIN	HEPARIN	SODIUM CITRATE 1:9	GEL	SODIUM CITRATE 1:4	EDTA	FLUORIDE

Weigao is committed to delivering a quality service and supply for blood collection products. Weigao works with the co-operation of hospitals and labs in order to ensure the correct result for the blood sample test; with Weigao precise and safe blood collection set, you can find the test result in a quick and precise way.

Heparin Tubes (Plasma)




Reference Number	Tube Size (mm)	Volume (ml)	Glass/ Plastic	Cap Colour	Pack/ Carton
SODIUM HEPARIN					
GS130752	13 x 75	2.0	G	GREEN	100 / 1800
PS120752	12 x 75	2.0	P	GREEN	100 / 1800
GS130753	13 x 75	3.0	G	GREEN	100 / 1800
PS120753	12 x 75	3.0	P	GREEN	100 / 1800
GS131005	13 x 100	5.0	G	GREEN	100 / 1200
PS121005	12 x 100	5.0	P	GREEN	100 / 1200
LITHIUM HEPARIN					
GL130752	13 x 75	2.0	G	GREEN	100 / 1800
PL120752	12 x 75	2.0	P	GREEN	100 / 1800
GL130753	13 x 75	3.0	G	GREEN	100 / 1800
PL120753	12 x 75	3.0	P	GREEN	100 / 1800
GL131005	13 x 100	5.0	G	GREEN	100 / 1200
PL121005	12 x 100	5.0	P	GREEN	100 / 1200

Plasma tube with Lithium Heparin, Sodium Heparin used in anti-coagulation not only for routine clinical biochemistry tests and emergency biochemistry tests but also for some test items in blood theology.

Reference centrifugation speed: 3000-3500 rpm for 5-8 mins

Additive status: spray



CE TÜV

EDTA Tubes (Whole Blood)




Reference Number	Tube Size (mm)	Volume (ml)	Glass/ Plastic	Cap Colour	Pack/ Carton
EDTA K2					
GE1307522	13 x 75	2.0	G	LAVENDER	100 / 1800
PE1207522	12 x 75	2.0	P	LAVENDER	100 / 1800
GE1307523	13 x 75	3.0	G	LAVENDER	100 / 1800
PE1207523	12 x 75	3.0	P	LAVENDER	100 / 1800
GE1307525	13 x 75	5.0	G	LAVENDER	100 / 1800
PE1210025	12 x 100	5.0	P	LAVENDER	100 / 1200
EDTA K3					
GE1307532	13 x 75	2.0	G	LAVENDER	100 / 1800
PE1207532	12 x 75	2.0	P	LAVENDER	100 / 1800
GE1307533	13 x 75	3.0	G	LAVENDER	100 / 1800
PE1207533	12 x 75	3.0	P	LAVENDER	100 / 1800
GE1307535	13 x 75	5.0	G	LAVENDER	100 / 1800
PE1210035	12 x 100	5.0	P	LAVENDER	100 / 1200

EDTA tubes can be used on whole blood hematology determinations, blood type verification and equipment.

1.8mg EDTA-K2 per 1ml volume

Additive status: Spray



CE TÜV

Plain Tubes (Serum)




Reference Number	Tube Size (mm)	Volume (ml)	Glass/Plastic	Cap Colour	Pack/ Carton
GP1307502	13 x 75	2.0	G	RED	100 / 1800
PP1207502	12 x 75	2.0	P	RED	100 / 1800
GP1307503	13 x 75	3.0	G	RED	100 / 1800
PP1207503	12 x 75	3.0	P	RED	100 / 1800
GP1307505	13 x 75	5.0	G	RED	100 / 1800
GP1310005	13 x 100	5.0	G	RED	100 / 1200
PP1210005	12 x 100	5.0	P	RED	100 / 1200
GP1310007	13 x 100	7.0	G	RED	100 / 1200
PP1210065	12 x 100	6.5	P	RED	100 / 1200
GP1610010	16 x 100	10.0	G	RED	100 / 1200
PP1610010	16 x 100	10.0	P	RED	100 / 1200

Serum tube can be used in collecting blood sample in clinical biochemistry and immunology.

Reference centrifugation speed: 3000-3500 rpm
Reference clotting time: No additive 40-60 mins,
Activator 10 mins

Additive status: Spray



GEL Tubes (Serum Separator)




Reference Number	Tube Size (mm)	Volume (ml)	Glass/Plastic	Cap Colour	Pack/ Carton
GG1307503	13 x 75	3.0	G	YELLOW	100 / 1800
PG1207503	12 x 75	3.0	P	YELLOW	100 / 1800
GG1310005	13 x 100	5.0	G	YELLOW	100 / 1200
PG1210005	12 x 100	5.0	P	YELLOW	100 / 1200
GG1610007	16 x 100	7.0	G	YELLOW	100 / 1200
PG1610007	16 x 100	7.0	P	YELLOW	100 / 1200
GG1610009	16 x 100	9.0	G	YELLOW	100 / 1200
PG1610009	16 x 100	9.0	P	YELLOW	100 / 1200

GEL can be used in separating serum and plasma, pre-filled for quick serum measurement and immunology test.

Reference centrifugation speed: 3000-3500rpm for 5-8MINS

Additive status: GEL




Sodium Citrate 1:4 Tubes (ESR)



Reference Number	Tube Size (mm)	Volume (ml)	Glass/ Plastic	Cap Colour	Pack/ Carton
GN1307516	13 x 75	1.6	G	BLACK	100 / 1800
PN1207516	12 x 75	1.6	P	BLACK	100 / 1800
GN1307524	13 x 75	2.4	G	BLACK	100 / 1800
PN1207524	12 x 75	2.4	P	BLACK	100 / 1800
GN0812016	8 x 120	1.6	G	BLACK	100 / 1200
GN0812012	8 x 120	1.28	G	BLACK	100 / 1200

Sodium Citrate 1:4 used with 3.2% (0.109 mol/L) sodium citrate for anti coagulation and ratio of anticoagulant and blood is 1:4

Ratio of anticoagulant and blood is 1:4
Additive status: liquid




Sodium Citrate 1:9 Tubes (Coagulation)



Reference Number	Tube Size (mm)	Volume (ml)	Glass/ Plastic	Cap Colour	Pack/ Carton
GN1307518	13 x 75	1.8	G	BLUE	100 / 1800
PN1207518	12 x 75	1.8	P	BLUE	100 / 1800
GN1307527	13 x 75	2.7	G	BLUE	100 / 1800
PN1207527	12 x 75	2.7	P	BLUE	100 / 1800
GN1307545	13 x 75	4.5	G	BLUE	100 / 1800
PN1207545	12 x 75	4.5	P	BLUE	100 / 1800

Sodium Citrate 1:9 with 3.8% (0.129mol/L) sodium citrate as anticoagulant used for coagulation test and ratio of anticoagulant and blood is 1:9.

Ratio of anticoagulant and blood is 1:9
Additive status: liquid



Fluoride Tubes (Blood Glucose)




Reference Number	Tube Size (mm)	Volume (ml)	Glass/ Plastic	Cap Colour	Pack/ Carton
GF1307502	13 x 75	2.0	G	GREY	100 / 1800
PF1207502	12 x 75	2.0	P	GREY	100 / 1800
GF1307503	13 x 75	3.0	G	GREY	100 / 1800
PF1207503	12 x 75	3.0	P	GREY	100 / 1800
GF1310005	13 x 100	5.0	G	GREY	100 / 1200
PF1210005	12 x 100	5.0	P	GREY	100 / 1200
GF1310007	13 x 100	7.0	G	GREY	100 / 1200
PF1210065	12 x 100	6.5	P	GREY	100 / 1200

Glucose Tube with additive of potassium oxalate/sodium fluoride EDTA is used for test of blood coagulation mechanism.


Reference centrifugation speed: 3000-3500 rpm for 5-8 mins

Additive status: Spray

CE TÜV



Blood Collection Needles




Reference Number	Needle Size	Length (mm)	Reference Colour	Type	Inner Box / Carton
VNP20100	20G	25	YELLOW	PEN LIKE	100 / 5000
VNP20112	20G	38	YELLOW	PEN LIKE	100 / 5000
VNB20100	20G	25	YELLOW	BUTTERFLY WINGS WITH TUBE	100 / 2000
VNP21100	21G	25	GREEN	PEN LIKE	100 / 5000
VNP21112	21G	38	GREEN	PEN LIKE	100 / 5000
VNB21100	21G	25	GREEN	BUTTERFLY WINGS WITH TUBE	100 / 2000
VNP22100	22G	25	BLACK	PEN LIKE	100 / 5000
VNP22112	22G	38	BLACK	PEN LIKE	100 / 5000
VNB22100	22G	25	BLACK	BUTTERFLY WINGS WITH TUBE	100 / 2000
VNH20000	PLASTIC NEEDLE HOLDER				100 / 2000

Vacuum blood collection needle and holder are used with vacuum tubes to collect the blood, the needle size is marked with different colors in 20G, 21G and 22G.

Note: Please use the needle by the professional user and under correct performance to avoid the needlestick injury.

CE 0123 TÜV

EQUINOX



EQUINOX X

Equinox Orthopaedics Pty Ltd
Level 35, 100 Barangaroo Ave,
International Tower 1, Barangaroo,
Sydney, NSW 2000
www.equinox-orthopaedics.com