

Solutions for Cannulated Screws



Ø3.0mm Cannulated Screw(Countersunk head/Short thread)




With drill bit: cannulated, ϕ 2.0mm
(WG REF #30402103)

Indication: It is used for the fracture on metacarpal, Carpus and Plantar phalanges.

Ref #	Size	Material
08.301.11.601133010	3.0x10mm	T
08.301.11.601133012	3.0x12mm	
08.301.11.601133014	3.0x14mm	
08.301.11.601133016	3.0x16mm	
08.301.11.601133018	3.0x18mm	
08.301.11.601133020	3.0x20mm	
08.301.11.601133022	3.0x22mm	
08.301.11.601133024	3.0x24mm	
08.301.11.601133026	3.0x26mm	
08.301.11.601133028	3.0x28mm	
08.301.11.601133030	3.0x30mm	
08.301.11.601133032	3.0x32mm	
08.301.11.601133034	3.0x34mm	
08.301.11.601133036	3.0x36mm	
08.301.11.601133038	3.0x38mm	
08.301.11.601133040	3.0x40mm	

Ø3.0mm Cannulated Screw(Countersunk head/half thread)



With drill bit: cannulated, ϕ 2.0mm
(WG REF #30402103)

Indication: It is used for the fracture on metacarpal, Carpus and Plantar phalanges.

Ref #	Size	Material
08.301.11.601153010	3.0x10mm	T
08.301.11.601153012	3.0x12mm	
08.301.11.601153014	3.0x14mm	
08.301.11.601153016	3.0x16mm	
08.301.11.601153018	3.0x18mm	
08.301.11.601153020	3.0x20mm	
08.301.11.601153022	3.0x22mm	
08.301.11.601153024	3.0x24mm	
08.301.11.601153026	3.0x26mm	
08.301.11.601153028	3.0x28mm	
08.301.11.601153030	3.0x30mm	
08.301.11.601153032	3.0x32mm	
08.301.11.601153034	3.0x34mm	
08.301.11.601153036	3.0x36mm	
08.301.11.601153038	3.0x38mm	
08.301.11.601153040	3.0x40mm	



Ø4.5mm Cannulated Screw (Countersunk head)



With drill bit: cannulated, Φ 3.2mm
(WG REF #30402700)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

Ref #	Size	Material
08.301.11.601354520	4.5x20mm	T
08.301.11.601354522	4.5x22mm	
08.301.11.601354524	4.5x24mm	
08.301.11.601354526	4.5x26mm	
08.301.11.601354528	4.5x28mm	
08.301.11.601354530	4.5x30mm	
08.301.11.601354532	4.5x32mm	
08.301.11.601354534	4.5x34mm	
08.301.11.601354536	4.5x36mm	
08.301.11.601354538	4.5x38mm	
08.301.11.601354540	4.5x40mm	
08.301.11.601354542	4.5x42mm	
08.301.11.601354544	4.5x44mm	
08.301.11.601354546	4.5x46mm	
08.301.11.601354548	4.5x48mm	
08.301.11.601354550	4.5x50mm	
08.301.11.601354552	4.5x52mm	
08.301.11.601354554	4.5x54mm	
08.301.11.601354556	4.5x56mm	
08.301.11.601354558	4.5x58mm	
08.301.11.601354560	4.5x60mm	
08.301.11.601354564	4.5x64mm	
08.301.11.601354568	4.5x68mm	
08.301.11.601354572	4.5x72mm	

Ø3.5mm Cannulated Screw (Ball head)



With drill bit: cannulated, ϕ 2.7mm
(WG REF #30400712)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

Ref #	Size	Material
08.301.11.600233510	3.5x10mm	T
08.301.11.600233512	3.5x12mm	
08.301.11.600233514	3.5x14mm	
08.301.11.600233516	3.5x16mm	
08.301.11.600233518	3.5x18mm	
08.301.11.600233520	3.5x20mm	
08.301.11.600233522	3.5x22mm	
08.301.11.600233524	3.5x24mm	
08.301.11.600233526	3.5x26mm	
08.301.11.600233528	3.5x28mm	
08.301.11.600233530	3.5x30mm	
08.301.11.600233532	3.5x32mm	
08.301.11.600233534	3.5x34mm	
08.301.11.600233536	3.5x36mm	
08.301.11.600233538	3.5x38mm	
08.301.11.600233540	3.5x40mm	
08.301.11.600233542	3.5x42mm	
08.301.11.600233544	3.5x44mm	
08.301.11.600233546	3.5x46mm	
08.301.11.600233548	3.5x48mm	
08.301.11.600233550	3.5x50mm	



Ø4.5mm Cannulated Screw (Ball head)




With drill bit: cannulated, ϕ 3.2mm
(WG REF #6090004)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

Ref #	Size	Material
08.301.01.600074020	4.5x20mm	T
08.301.01.600074022	4.5x22mm	
08.301.01.600074024	4.5x24mm	
08.301.01.600074026	4.5x26mm	
08.301.01.600074028	4.5x28mm	
08.301.01.600074030	4.5x30mm	
08.301.01.600074032	4.5x32mm	
08.301.01.600074034	4.5x34mm	
08.301.01.600074036	4.5x36mm	
08.301.01.600074038	4.5x38mm	
08.301.01.600074040	4.5x40mm	
08.301.01.600074042	4.5x42mm	
08.301.01.600074044	4.5x44mm	
08.301.01.600074046	4.5x46mm	
08.301.01.600074048	4.5x48mm	
08.301.01.600074050	4.5x50mm	
08.301.01.600074052	4.5x52mm	
08.301.01.600074054	4.5x54mm	
08.301.01.600074056	4.5x56mm	
08.301.01.600074058	4.5x58mm	
08.301.01.600074060	4.5x60mm	
08.301.01.600074064	4.5x64mm	
08.301.01.600074068	4.5x68mm	
08.301.01.600074072	4.5x72mm	

Ø7.3mm Cannulated Screw (Ball head/long thread)



With drill bit: cannulated, ϕ 5.0mm
(WG REF #6090202)

Indication: It is used for the fracture on femur neck and the fusion on tibia and astragalus and calcaneus.

Ref #	l	Size	Material
08.301.01.600058050	32	7.3x50mm	T
08.301.01.600058055	32	7.3x55mm	
08.301.01.600058060	32	7.3x60mm	
08.301.01.600058065	32	7.3x65mm	
08.301.01.600058070	32	7.3x70mm	
08.301.01.600058075	32	7.3x75mm	
08.301.01.600058080	32	7.3x80mm	
08.301.01.600058085	32	7.3x85mm	
08.301.01.600058090	32	7.3x90mm	
08.301.01.600058095	32	7.3x95mm	
08.301.01.600058100	32	7.3x100mm	
08.301.01.600058105	32	7.3x105mm	
08.301.01.600058110	32	7.3x110mm	
08.301.01.600058115	32	7.3x115mm	
08.301.01.600058120	32	7.3x120mm	

Ø7.3mm Cannulated Screw (Ball head/short thread)




With drill bit: cannulated, ϕ 5.0mm
(WG REF #6090202)

Indication: It is used for the fracture on femur neck and the fusion on tibia and astragalus and calcaneus.

Ref #	l	Size	Material
08.301.01.600048050	16	7.3x50mm	T
08.301.01.600048055	16	7.3x55mm	
08.301.01.600048060	16	7.3x60mm	
08.301.01.600048065	16	7.3x65mm	
08.301.01.600048070	16	7.3x70mm	
08.301.01.600048075	16	7.3x75mm	
08.301.01.600048080	16	7.3x80mm	
08.301.01.600048085	16	7.3x85mm	
08.301.01.600048090	16	7.3x90mm	
08.301.01.600048095	16	7.3x95mm	
08.301.01.600048100	16	7.3x100mm	
08.301.01.600048105	16	7.3x105mm	
08.301.01.600048110	16	7.3x110mm	
08.301.01.600048115	16	7.3x115mm	
08.301.01.600048120	16	7.3x120mm	



Ø3.0mm Self- tapping Cannulated Screw



With drill bit: Φ 2.0 Cannulated drill (WG REF #30402206)
 Φ 2.0 conical drill (WG REF# 30402506)

Indication: It is used for the fracture on metacarpal ,
 Carpus and Plantar phalanges.

Ref #	Size	Material
08.301.11.601573008	3.0x8mm	T
08.301.11.601573010	3.0x10mm	
08.301.11.601573012	3.0x12mm	
08.301.11.601573014	3.0x14mm	
08.301.11.601573016	3.0x16mm	
08.301.11.601573018	3.0x18mm	
08.301.11.601573020	3.0x20mm	

Ø3.5mm Self- tapping Cannulated Screw




With drill bit: Φ 2.5 cannulated bone drill (WG REF #30402106)
 Φ 2.5 Conical bone drill (WG REF #30402406)

Indication: It is used for the fracture on metacarpal ,
 Carpus and Plantar phalanges.

Ref #	Size	Material
08.301.11.601613516	3.5x16mm	T
08.301.11.601613518	3.5x18mm	
08.301.11.601613520	3.5x20mm	
08.301.11.601613522	3.5x22mm	
08.301.11.601613524	3.5x24mm	
08.301.11.601613526	3.5x26mm	
08.301.11.601613528	3.5x28mm	
08.301.11.601613530	3.5x30mm	

Ø4.0mm Self- tapping Cannulated Screw




With drill bit: Φ 3.0 canulated bone drill (WG REF #30402006)
 Φ 3.0 Conical bone drill (WG REF #30402306)

Indication: It is used for the fracture on metacarpal , Carpus and Plantar phalanges.

Ref #	Size	Material
08.301.11.601654016	4.0x16mm	T
08.301.11.601654018	4.0x18mm	
08.301.11.601654020	4.0x20mm	
08.301.11.601654022	4.0x22mm	
08.301.11.601654024	4.0x24mm	
08.301.11.601654026	4.0x26mm	
08.301.11.601654028	4.0x28mm	
08.301.11.601654030	4.0x30mm	

Ø4.5mm Self- tapping Cannulated Screw




With drill bit: Φ 4.6- 3.2 drill (WG REF #30402906)
 Φ 4.4- 3 drill (WG REF #30403006)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

Ref #	Size	Material
08.301.11.601694524	4.5x24mm	T
08.301.11.601694526	4.5x26mm	
08.301.11.601694528	4.5x28mm	
08.301.11.601694530	4.5x30mm	
08.301.11.601694532	4.5x32mm	
08.301.11.601694534	4.5x34mm	
08.301.11.601694536	4.5x36mm	
08.301.11.601694538	4.5x38mm	
08.301.11.601694540	4.5x40mm	
08.301.11.601694542	4.5x42mm	
08.301.11.601694544	4.5x44mm	
08.301.11.601694546	4.5x46mm	
08.301.11.601694548	4.5x48mm	
08.301.11.601694550	4.5x50mm	



Ø5.0mm Self- tapping Cannulated Screw




With drill bit: Φ 4-5.6 drill (WG REF #30403606)
 Φ 4.5- 6 drill (WG REF # 30403706)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

Ref #	Size	Material
08.301.11.601735035	5.0x35mm	T
08.301.11.601735040	5.0x40mm	
08.301.11.601735045	5.0x45mm	
08.301.11.601735050	5.0x50mm	
08.301.11.601735055	5.0x55mm	
08.301.11.601735060	5.0x60mm	
08.301.11.601735065	5.0x65mm	
08.301.11.601735070	5.0x70mm	
08.301.11.601735075	5.0x75mm	
08.301.11.601735080	5.0x80mm	

Ø7.0mm Self-tapping Cannulated Screw



With drill bit: Φ 4.5- 6.7
 Reamer 1 (WG REF#30404406)
 Reamer 2 (WG REF#30404506)
 Reamer 3 (WG REF#30404606)
 Φ 5.0- 6.7
 Reamer 1 (WG REF# 30404706)
 Reamer 2 (WG REF# 30404806)
 Reamer 3 (WG REF# 30404906)

Indication: It is used for the fusion of tibia and astragalus.

Ref #	Size	Material
08.301.11.601777040	7.0x40mm	T
08.301.11.601777045	7.0x45mm	
08.301.11.601777050	7.0x50mm	
08.301.11.601777055	7.0x55mm	
08.301.11.601777060	7.0x60mm	
08.301.11.601777065	7.0x65mm	
08.301.11.601777070	7.0x70mm	
08.301.11.601777075	7.0x75mm	
08.301.11.601777080	7.0x80mm	
08.301.11.601777085	7.0x85mm	
08.301.11.601777090	7.0x90mm	
08.301.11.601777095	7.0x95mm	
08.301.11.601777100	7.0x100mm	