






# Broad Range of Extraction Kits for Respiratory, Oncology & Fertility

<b>Extraction Kit</b>	<b>32A/B Viral RNA Extraction Kit</b>	<b>Pathogenic RNA/DNA Extraction Kit</b>	<b>ReadiMag FFPE DNA or RNA Extraction Kit</b>	<b>ReadiMag cfDNA Extraction Kit</b>	<b>ReadiMag gDNA Extraction Kit</b>
<b>System Compatibility</b>	ANDiS 350/360	ANDiS 350/360	ANDiS 350/360	ANDiS 360	ANDiS 360
					
<b>Sample type</b>	Nasopharyngeal, Oropharyngeal, Saliva, Cervical	Nasopharyngeal, Oropharyngeal, Saliva, Sputum, Cervical & Urethral Swab, Urine, Stool	FFPE tissue	Serum, plasma	Saliva, whole blood and buffy coat
<b>Run time</b>	<15 minutes (Dependant on sample types and instrument used)	<15 minutes (Dependant on sample types and instrument used)	45 minutes	45 minutes	50 minutes
<b>Sample input volume</b>	200-400µL	<400µL	5-40µm	2000-4000µL	<200µL
<b>Elution volume</b>	100µL	100µL	50µL	50µL	50µL

