


cfDNA Extraction Kit for NIPT and other Fertility Screenings

- 3DMed's cfDNA extraction kit features pre-filled reagent plates, that simplify the extraction protocol by reducing the need for manual handling
- Specialised plates allow for the extraction of high-quality cfDNA from high-volume, low-concentration samples
- Application for extraction of cfDNA in NIPT



ReadiMag cfDNA Extraction Kit	
Pack size	16 tests/kit
Sample type	Serum, plasma
Processing time	45 minutes (+10 min pre-processing)
Sample volume	2000-4000 μ L
Elution volume	50 μ L
System compatibility	ANDiS 360
Areas of application	NIPT, NGS, liquid biopsy