


## gDNA Extraction Kit for NGS including Oncology

- 3DMed's gDNA extraction kit features pre-filled reagent plates, that simplify the extraction protocol by reducing the need for manual handling
- Application for the potential extraction of gDNA for preimplantation screening and reproductive carrier screening



ReadIMag gDNA Extraction Kit	
Pack size	64 tests/kit
Sample type	Blood, saliva
Processing time	50 minutes
Sample volume	200µL
Elution volume	50µL
System compatibility	ANDiS 360
Areas of application	NGS, Fertility screening (validation required)