




Extraction Kit	cfDNA Extraction Kit	ReadiMag FFPE DNA Extraction Kit	ReadiMag FFPE RNA Extraction Kit
System Compatibility	<i>ANDiS 360</i>	<i>ANDiS 350/360</i>	<i>ANDiS 350/360</i>
			
Sample type	Serum, plasma	FFPE tissue	FFPE tissue
Run time	45 minutes	45 minutes	45 minutes
Sample input volume	2000-4000 μ L	5-40 μ m	5-40 μ m
Elution volume	50 μ L	50 μ L	50 μ L

